

# SEDEX LT-ELSD™

## THE RESULT OF 20 YEARS OF EVOLUTION



**SEDEX<sup>LT</sup>-ELSD™**  
Low Temperature  
Evaporative Light Scattering Detectors



SENSITIVITY  
FLEXIBILITY  
EXPERIENCE

**SENSITIVITY**  
**FLEXIBILITY**  
**EXPERIENCE**



## A BETTER WAY TO GET MORE FROM HPLC AND LC/MS ANALYSIS

### Introduction to SEDEX DETECTORS

SEDERE develops, manufactures, distributes and supports SEDEX detectors, the most complete product line dedicated to LT-ELSD™. As one of the originators of LT-ELSD, SEDERE remains exclusively focused on this technology as a core competency.

As the industry leader, SEDERE leverages decades of experience and customer knowledge to continually raise the bar for sensitive detector performance in chromatography laboratories.

To keep up with evolving user needs, SEDERE has added two new detectors, the SEDEX Models 80LT and 85LT. The unparalleled selection of LT-ELSD detection systems can match performance requirements and budgets for all chromatography applications.

### The challenges for sensitive, complete LT-ELSD™ require specialized design and knowledge of HPLC APPLICATIONS

#### FEATURES

- Applicable to semi-volatile and thermo-sensitive compounds
- Lowest background noise to provide excellent S/N ratio
- Optimization of peak shape and peak width
- Consistency of operating protocols
- Compatibility of nebulization with any HPLC protocol
- Prevent contamination of critical detector components
- User friendly, low maintenance system
- Integrate readily with HPLC software

#### SEDEX TECHNOLOGY

- The strength of the real Low Temperature technology
- An enhanced digital signal processing
- Nebulizer design and data rate up to 100Hz
- Complete, efficient and reliable information and SOPs
- Four nebulizers to cover a complete range of flow-rates
- Safety features, patented Gas Supported Focusing (GSF™)
- Plug-and-play detector; power-down methods
- RS 232

*Only SEDEX meets all the challenges.*

## POWER HPLC WITH ENRICHED DETECTION


### The Reasons to Adopt LT-ELSD™

- Provides an essentially universal measurement **not** dependent on a compound's:
  - Absorbance
  - Electroactivity
  - Fluorescence
  - Radioactivity etc.
- Provides mass detector capability
- Minimizes volatilization of compounds through Low Temperature evaporation
- Allows gradient elution
- Permits a wide choice of volatile buffers
- Works with underivatized compounds
- Surpasses RI sensitivity by orders of magnitude

### Specifically, it is a Technology Solution

for detecting many compounds critical to pharmaceutical products!:

- Phospholipids, amino acids
- Polysorbates
- Polyethylene glycols (PEGs)
- Stearic, oleic, citric, tartaric acids
- Corn oil, castor oil, coconut oil
- Carbohydrates






*Carrageenan Analysis  
by LT-ELSD™*

**LT-ELSD™ PREVENTS YOU FROM MISSING  
COMPOUNDS "INVISIBLE" TO UV/VIS DETECTORS.**

# SENSITIVITY

SEDEX<sup>LT</sup>-ELSD™ Low Temperature Evaporative Light Scattering Detectors

## THREE STAGES OF SEDEX LT-ELSD EACH OPTIMIZED FOR HIGH PERFORMANCE DETECTION



### Nebulize Eluent and Select Small Droplets to minimize background noise


The eluent from the HPLC column is transformed into a fine mist of droplets, using the nebulizer. The SEDEX LT technology allows the selection of droplets as a function of their size in order to prevent larger droplets from entering the evaporation (drift) tube. Large droplets are responsible for increased background noise, as they are more difficult to evaporate. This selection of droplets enables detection using a very low evaporation temperature, with resulting low baseline noise and excellent sensitivity to solutes, even those which are semi-volatile.

### Evaporate at Low Temperature every time so you won't miss any compound


Solute molecules are obtained from the mist using a heated evaporation (drift) tube, using the Low Temperature technology. All SEDEX detectors are designed to evaporate high boiling mobile phases (e.g. those with high aqueous concentration) at very low temperatures. This unique low temperature provided in the standard operating mode minimizes the potential for thermal decomposition of the compound of interest, and makes the SEDEX method a more reliable way to detect everything in the sample.

### Detect Light Scattering using Patented Gas Supported Focusing (GSF™) for less maintenance and better data

The solute molecules from the mist, assisted by Gas Supported Focusing (GSF), pass through an optical head designed to measure light scattering. Gas supported focusing (GSF) involves the addition of gas to focus the solute particles within the optical head for enhanced detection. The magnitude of resulting scattered light is related to the concentration of the compounds in the sample (i.e. mass). The optical head is also protected from contamination by the gas surrounding the solute particles.



- SEDEX LT-ELSD™ eliminates large droplets to significantly lower baseline noise.
- The quantity of droplets evaporated is greater with LT-ELSD than with other detectors using splitter methods, providing enhanced sensitivity.



The SEDEX Low Temperature technology is the key to significantly better sensitivity than other ELSD systems that may use splitters or require high temperature to reduce noise.

SENSITIVITY

FLEXIBILITY

EXPERIENCE

## SEDEX: MORE QUALITY, MORE OPTIONS

### Introducing the Visually Enhanced series with SEDEX 80 LT-ELSD™, the New Standard

SEDEX 80 LT-ELSD™ combines high sensitivity, reliability and accuracy for your analyses, using the thanks to unrivalled SEDEX Low Temperature technology. This detector presents a new **Visually Enhanced** design and a number of outstanding innovations providing the best optical and electronic benefits at a **very competitive price**.

Two different nebulizers are available to optimize sensitivity and resolution at HPLC flow rates from 100µL/min to 5mL/min.

The SEDEX 80 LT-ELSD™ is user-friendly and meets the requirements of any HPLC system and application. In addition, the detector features programmable power-down methods and optimized controls.



SEDEX 80LT




SEDEX 75LT

### SEDEX 75 LT-ELSD™, the Reference Detector

- More than 2000 Sedex 75 are employed in the field today.
- Includes Low Temperature technology and patented GSFT™ (Gas Supported Focusing).
- Four specific nebulizers are available to optimize sensitivity and resolution at HPLC flow rates from 5µL/min to 5mL/min.
- Easy to use and compatible with all HPLC systems.

## SEDEX NEBULIZER OPTIONS





# FLEXIBILITY

## Introducing the Visually Enhanced series with SEDEX 85 LT-ELSD™, the State-of-the-Art

The new **Visually Enhanced** SEDEX 85 LT-ELSD™ employs the most advanced technology to provide the highest sensitivity, efficiency and reproducibility of all ELS detectors in the market. This powerful and versatile instrument can be used with any conventional HPLC, from micro to preparative, due to the use of four specific nebulisers. In addition, it is optimized to meet the requirements of ultra fast liquid chromatography including Ultra High Performance Liquid Chromatography (U-HPLC), High Temperature Liquid Chromatography (HTLC) and  $\mu$ -HTLC. For  $\mu$ -HTLC, a patented cell has been designed to provide the lowest dispersion and the best peak shape. As an example of the power of the system, the SEDEX 85 can measure **sub-one second** peak widths with data rate up to 100Hz.

As a true mass detector, calibration using SEDEX 85 LT-ELSD™ is straightforward and consistent in responses with variations below 10% between compounds having identical concentrations and belonging to the same chemical class. In addition, linearity over 3 orders of magnitude can be easily achieved using the log-log coordinates of the general equation  $A=a.m^b$ . Other equations such as second-order or third-order polynomial equations can also provide excellent quantitative results with correlation coefficients nearly equal to 1.



SEDEX 85LT

## THE OPTIMIZED PERFORMANCE OF SEDEX 85


### FEATURES

- Improved optical head design
- Digital signal processing
- Data rate 100Hz
- Software controlled
- Option of four interchangeable nebulizers
- Visual enhancement


### BENEFITS

- Lowest limits of detection: 100pg level sensitivity
- Contamination free
- Minimizes peak broadening, improves Signal/Noise ratio
- Allows **sub-one second** peak width
- Provides full remote control including gain programming
- Integrable in any HPLC software method, GLP compliant
- Assures synergy with HPLC methods improves sensitivity, peak efficiency and peak symmetry
- Controls visually nebulization

### U-HPLC: Paraben Analysis



### HTLC: Acetamidophenol, Hydroxydibenzofuran, tert-Butylphenoxybenzaldehyde




## CASE STUDIES:

### CASE 1 PHOSPHOLIPIDS<sup>2</sup>

LT-ELSD™ solves the major problems common to other HPLC detectors: lack of sensitivity, incompatibility with multi-solvent gradients... This state-of-the-art technique is ideally suited to non-chromophoric compounds, such as lipids and phospholipids, which have a low extinction coefficient.

LT-ELSD™ is also highly useful where the mobile phase contains a chromophore, which blanks out the UV detector.


**Column:** Lichrospher™ Si 60, 5µm  
**Injected Sample:** 3µg quantity, Phospholipids  
**Eluent:** CHCl<sub>3</sub>/MeOH/H<sub>2</sub>O/30%NH<sub>3</sub>  
**Gradient:** 60/34/5.5/0.5 to 80/19.5/0.5/0.0 in 14 min then hold for 10 min  
**Flow Rate:** 1mL/min



- |                                      |
|--------------------------------------|
| 1. Apolar lipids                     |
| 2. CR Cerebroside - Not hydroxylated |
| 3. CROH Cerebroside - Hydroxylated   |
| 4. PE Phosphatidyl ethanolamine      |
| 5. PI Phosphatidyl inositol          |
| 6. PS Phosphatidyl serine            |
| 7. PC Phosphatidyl choline           |
| 8. PA Phosphatidic acid              |
| 9. Sph Sphingomyeline                |
| 10. LPC Lysophosphatidyl choline     |

### CASE 2 MONO AND OLIGOSACCHARIDES

Unlike RI Detection, LT-ELSD™ allows gradient elution. Gradient elution provides increased resolution of sugars in minimal time, impossible with RI and isocratic elution. Moreover, lower detectable limits (sensitivity) can improve by orders of magnitude. Nanomole and picomole detectability are obtained with the improved sensitivity of LT-ELSD™. In spite of a complex matrix, these mono and oligosaccharides are easily and rapidly characterized by gradient HPLC with LT-ELSD™. Previously, RI detection entailed slow and tedious programmed flow, often up to several hours. LT-ELSD™ also enables analysis of high "DPs" which is an important advantage.




**Column:** Premier Carbohydrate, 5µm, 250x4.6mm  
**Injected Sample:** Extract from minced onions  
**Eluent:** - A: MeCN  
          - B: Water + 0.0004N NH<sub>4</sub>OH  
**Gradient:** From 60 to 65% B in 15 min  
**Flow Rate:** 1mL/min

- |             |                              |
|-------------|------------------------------|
| 1. Fructose | 6. Nystose                   |
| 2. Glucose  | 7. DP4                       |
| 3. Sucrose  | 8. Fructofuranosyl-D-nystose |
| 4. DP3      | 9. DP5                       |
| 5. Kestose  |                              |

### CASE 3 UNDERIVATIZED AMINO ACIDS<sup>3</sup>


Analysis of amino acids has typically been complicated by the absence of adequate visible or ultraviolet chromophores in naturally occurring amino acids. Using LT-ELSD™, sensitivity is excellent with detection limits as low as 200 picomoles. In this study, twenty protein amino acids have been separated and quantified within 40min without any sample preparation step for derivatization.



**Column:** Hypercarb™ S, 5µm, 100x2.1mm  
**Injected Sample:** Mixture of 20 underivatized amino acids  
**Eluent:** - A: Water + NFPA 20mM  
          - B: MeCN  
**Gradient:** 0 to 15% B in 10 min, 15 to 26% B in 10 min, 26 to 50% B in 10 min  
**Flow Rate:** 0.2mL/min

### CASE 4 UNDERIVATIZED PEPTIDES<sup>4</sup>

In peptide "mapping", where gradient elution is required, LT-ELSD™ has a key advantage over UV detection: its baseline is unperturbed by the mobile phase change during the gradient, and remains flat. As a mass detector, LT-ELSD™ can also provide a material balance purity analysis. Degradation products are often lacking of the chromophores that were initially present in the original compounds, and do not allow purity evaluation using a UV detector:




**Column:** Hypercarb™ S, 5µm, 100x2.1mm  
**Injected Sample:** Mixture of 13 peptides  
**Eluent:** - A: Water + NFPA 20mM  
          - B: MeCN  
**Gradient:** 0 to 50% B in 60 min  
**Flow Rate:** 0.2mL/min

## CASE 5

### NATURAL PRODUCTS<sup>5</sup>

Many natural products such as herbal drugs are gaining more and more interest in the pharmaceutical and nutraceutical industry because they contain bioactive compounds. Some of these compounds such as saponins and terpenes do not possess any chromophore and therefore cannot be analyzed in HPLC using a UV detector. LT-ELSD™ can detect chromophoric and non-chromophoric molecules in a single gradient HPLC analysis with an excellent sensitivity, thanks to the true Low Temperature technology. This example shows a validated method for a simultaneous determination of isoflavones and saponins in soybean.




**Column:** Hypersil C18, 5 μm, 150x4.6mm  
**Injected Sample:** Soybean extract  
**Eluent:** - A: H<sub>2</sub>O + 0.025% TFA  
 - B: MeCN + 0.025% TFA  
**Gradient:** 85:15 to 60:40 in 15 min  
 and 60:40 to 45:55 in 20 min  
**Flow Rate:** 1.0mL/min

1 to 6. Isoflavones
7 to 10. Saponins

## CASE 6

### POLYETHYLENE GLYCOL AND TRITON

The high sensitivity and time saving potential of LT-ELSD™ are evident in the HPLC/LT-ELSD™ analysis of two polymer mixtures in a single run; this is not feasible with competing detections such as RI, UV and MS.




**Column:** Hypercarb™ S, 7 μm, 100x4.6mm  
**Injected Sample:** Mixture containing PEG400 and Triton X100  
**Eluent:** H<sub>2</sub>O / MeCN / CH<sub>2</sub>Cl<sub>2</sub>  
**Gradient:** From 80:20:0 to 0:100:0 in 15 min  
 From 0:100:0 to 0:20:80 in 25 min  
**Flow Rate:** 1mL/min

## CASE 7

### INORGANIC IONS

LT-ELSD™ can dramatically simplify the analysis of ions in aqueous samples. A broad range of volatile buffers can be used to separate the ions<sup>6</sup>. Since the mobile phase and buffers are vaporized before the ions are detected, the need for ion suppression is eliminated. This example shows an outstanding patented method<sup>7</sup> to determine immediately and simultaneously cations and anions in a mineral water.




**Column:** Hypercarb™ S, 7 μm, 100x4.6mm,  
 + Lichrosil IC CA, 5 μm, 100x4.6mm  
**Injected Sample:** Mineral water  
**Eluent:** - A: H<sub>2</sub>O  
 - B: HCOOH 100mM / NH<sub>4</sub>OH 60mM, pH : 3.71  
**Gradient:** 25%B during 7 min,  
 from 25 to 100%B in 8 min  
**Flow Rate:** 0.8mL/min

1. unknown	5. K <sup>+</sup>
2. Cl <sup>-</sup>	6. NO <sub>3</sub> <sup>-</sup>
3. Na <sup>+</sup>	7. Mg <sup>2+</sup>
4. SO <sub>4</sub> <sup>2-</sup>	8. Ca <sup>2+</sup>

## CASE 8

### U-HPLC/ELSD<sup>8</sup>

The pharmaceutical discovery environment requires an increasing number of rapid high throughput methods such as U-HPLC to determine the identity, purity and quantity of small molecules. In this regard, the powerful and versatile LT-ELSD™ is the detector of choice because of its universality, highest sensitivity and optimized technology which provides the smallest peak widths, the best symmetry and the highest data rate. This example shows an application which combines an ultra fast liquid chromatography system with LT-ELSD™, to determine the non-chromophoric artemisinin and derivatives used as antimalaric drugs.



**Column:** Acquity C18, 1.7 μm, 30x2.1mm  
**Injected Sample:** Artemisinin and derivatives  
**Eluent:** H<sub>2</sub>O / MeCN / TFA 38:60:2  
**Flow Rate:** 0.4mL/min

1. unknown
2 to 7. Artemisinin (6) and derivatives

SPECIFICATIONS	SEDEX 75LT	SEDEX 80LT	SEDEX 85LT
<b>COMPONENTS</b>			
Detection	Photomultiplier (PMT)	Photomultiplier (PMT)	Photomultiplier (PMT)
Light Source	Halogen Lamp	LED (470nm) Elapsed Time Counter	LED (470nm) Elapsed Time Counter Normalization
Temperature Range	Ambient - 100°C	Ambient - 100°C	Ambient - 100°C
Nebulizer	4 options: Capillary, Low Flow, HPLC, CC	2 options: HPLC, CC	4 options: Capillary, Low Flow, HPLC, CC
Eluent Flow Rate	5µL/min - 5mL/min	100µL/min - 5mL/min	5µL/min - 5mL/min
<b>DATA</b>			
Analog Output	0 - 1V	0 - 1V	0 - 1V
Gain Settings	2 <sup>11</sup> (2048)	2 <sup>11</sup> (2048)	2 <sup>11</sup> (2048)
Filter	Time Constant (None, Medium, High)	Moving Average (0 - 10s)	Moving Average (0 - 10s)
Data Rate	N/A	40Hz	100Hz
<b>COMMUNICATION</b>			
Selection & Display	Liquid Cristal Display and Keypad Windows based PC control	Liquid Cristal Display and Keypad	Liquid Cristal Display and Keypad Windows based PC control
Event	Contact closure for auto-zero	Contact closure for ready, auto-zero	Contact closure for ready, auto-zero
Power-down Methods	N/A	Shut-off: Gas, LED, Heating and/or PMT Cleaning Mode	Shut-off: Gas, LED, Heating and/or PMT Cleaning Mode
Computer Interface	RS232	RS232	RS232
Software	Option	Option	Bundled - Control of all the parameters (Windows 9x, NT, 2000, XP)
GLP Compliance	-	-	Audit trail - Password - Admin. mode
<b>EXTERNAL REQUIREMENTS</b>			
Power	115V/60Hz, 1.5A 230V/50Hz, 1.7A	115V/60Hz, 1.8A 230V/50Hz, 1.7A	115V/60Hz, 1.8A 230V/50Hz, 1.7A
Gas Pressure (Nitrogen or Air)	3.5 bars (51psi)	3.5 bars (51psi)	3.5 bars (51psi)
Dimensions	360mm (14.1") W 504mm (19.8") H 500mm (19.7") D	250mm (10") W 480mm (18") H 550mm (22") D	250mm (10") W 480mm (18") H 550mm (22") D
Weight	20Kg (44lbs)	18.5Kg (40lbs)	18.5Kg (40lbs)

## SEDERE is committed to user satisfaction with every SEDEX Detector, and supports you with :

- On-site installation and user training
- Full SOPs (Standard Operating Procedures) including IQ, OQ, PQ
- Formal Operational Qualification
- Technical and applications support
- Web-access to applications in many fields
- User seminars, on and off-site
- Flexible service contract options
- Easy-to-order spare parts and accessories
- Nebulizer options for a wide range of applications
- Web-based software upgrades



### SEDERE S.A.

BP 27 - Parc Volta - 94141 Alfortville Cedex FRANCE  
Tel: +33 (0)1 45 18 05 18 Fax: +33 (0)1 45 18 05 25  
Email: manager@sedere.com

[www.sedere.com](http://www.sedere.com)

### SEDERE INC.

4 Princess Road – Suite 201  
Lawrenceville, NJ 08648 USA  
Tel : (888) 452-4253 Fax : (607) 895-7158



# EXPERIENCE

"ArQule has been using SEDEX LT-ELSD™ for over five years to fully characterize automated parallel synthesis products. We have worked closely with the manufacturer, SEDERE, to adapt the instruments to the changing needs of high-speed analysis. In our laboratory, each SEDEX detector analyzes tens of thousands of samples every year in a 24/7 operation.

The instruments have proven to be robust and reliable in our high throughput environment."

– Wolfgang K. Goetzinger, Ph.D.  
Director of Analytical Chemistry  
ArQule, Inc.

## Ordering information

<b>SEDEX 85 LT-ELSD™</b> .....	Part #
HPLC Version 230 V .....	85000
HPLC Version 115V .....	85001
LF Version 230V .....	85300
LF Version 115V .....	85301
CC Version 230 V .....	85400
CC Version 115 V .....	85401
Capillary Version 230 V .....	85600
Capillary Version 115 V .....	85601

<b>SEDEX 80 LT-ELSD™</b> .....	Part #
HPLC Version 230V .....	80000
HPLC Version 115V .....	80001
CC Version 230V .....	80400
CC Version 115V .....	80401

<b>ACCESSORIES</b> .....	Part #
Gas regulator with 0.01µm filter and manometer .....	45100
TTL to CC converter .....	85080
Nitrogen generator .....	30100

<b>SEDEX 75 LT-ELSD™</b> .....	Part #
HPLC Version 230 V .....	75000
HPLC Version 115V .....	75001
LF Version 230V .....	75300
LF Version 115V .....	75301
CC Version 230 V .....	75400
CC Version 115 V .....	75401
Capillary Version 230 V .....	75600
Capillary Version 115 V .....	75601
HPLC/SFC Version 230 V .....	75200
HPLC/SFC Version 115 V .....	75201

<b>SEDEX 754 LT-ELSD</b> .....	Part #
Multichannel Version 230 V .....	754000-220
Multichannel Version 115 V .....	754000-115

# EXPERIENCE

## SEDEX<sup>LT</sup>-ELSD<sup>TM</sup>

### An Industry Standard for Low Temperature Evaporative Light Scattering Detection

The arrival of the Ultra Fast HPLC has fueled the demand for technology capable of both qualitative and quantitative analysis of complex mixtures at high speed. SEDEX LT-ELSD<sup>TM</sup> technology has been validated by extensive applications within the drug discovery, pharmaceutical and nutraceutical industries. SEDEX detectors are used in every major pharmaceutical company and in hundreds of biotechnology laboratories in industry and universities.

For many research and process requirements, complementary detection by SEDEX LT-ELSD<sup>TM</sup> has proven indispensable to high quality LC/MS and other HPLC procedures. SEDEX LT-ELSD<sup>TM</sup> is particularly valuable for effective compound library screening, where sample characterization may be incomplete. With other ELSD detectors, volatilization could limit the detection capability of the platform, resulting in loss of vital data.

By combining reliability and sensitivity, SEDEX detectors have taken their place in the armamentarium of excellent techniques for medicinal chemistry.

- 1- Lafosse, M. ; Elfakir, C. ; Morin-Allory, L. ; Dreux, M. The advantages of evaporative light scattering detection in pharmaceutical analysis by high performance liquid chromatography and supercritical fluid chromatography. *J. High Resolut. Chromatogr.* 1992, 15, 312.
- 2- Becart, J. ; Chevalier, C. ; Biesse, J.P. Quantitative analysis of phospholipids by HPLC with a light-scattering evaporating detector: Application to raw materials of cosmetic use. *J. High Resol. Chromatogr.*, 1990, 13, 126.
- 3- Chaimbault, P. ; Petritis, K. ; Elfakir, C. ; Dreux, M. Ion-pair chromatography on a porous graphitic carbon stationary phase for the analysis of twenty underivatized protein amino acids. *J. Chromatogr. A*, 2000, 870, 245.
- 4- Adoubel, A.A. ; Guenu, S. ; Elfakir, C. ; Dreux, M. Separation of underivatized small peptides on a porous graphitic carbon column by ion-pair chromatography and evaporative light-scattering detection. *J. Liq. Chrom. & Rel. Technol.*, 2000, 16, 2433.
- 5- Ganzera, M. ; Stuppner, H. ; Khan, I. Simultaneous Determination of Saponins and Isoflavones in Soybean (*Glycine max L.*) by Reverse-Phase Liquid Chromatography with Evaporative Light-Scattering and Ultraviolet Detection. *J. AOAC Intern.*, 2004, 87, 1189.
- 6- Petritis, K. ; Dessans, H. ; Elfakir, C. ; Dreux, M. Volatility evaluation of mobile-phase/ electrolyte additives for mass spectrometry. *LC GC Eur.* 2001, 15, 98.
- 7- United States Patent - Number: 6,148,661 - Date: 11/21/2000 - Method of Separating and Rapidly Analyzing a Sample - Inventor: Michel Dreux, SEDERE SA Alfortville and Université d'Orléans, France.
- 8- Russo, R. ; Guillaume, D. ; Bicchi, C. ; Rudaz S. ; Veuthey, J.L. UPLC coupled to an evaporative light scattering detector for the rapid analysis of non UV detectable compounds. Poster for SEP07 (2007, 20-22 March), Grenoble, France.